aura B

aura B 24

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ITALIANO

Instructions for use

ENGLISH

DXYA106 Rev_6

Costruito da / The Manufacturer

MEDILINE ITALIA s.r.l.

Via 8 marzo, 4
43025 Corte Tegge – Cavriago (RE)
ITALY

Gentile Cliente,

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Product Page: DAP-1160
Firmware Version: v1.20

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Setup Wizard
 Wireless Setup
LAN Setup
 Logout

NETWORK SETTINGS :

Use this section to configure the internal network settings of your AP and also to configure the built-in DHCP Server to assign IP addresses to the computers on your network. The IP Address that is configured here is the IP Address that you use to access the Web-based management interface. If you change the IP Address here, you may need to adjust your PC's network settings to access the network again.

LAN CONNECTION TYPE :

Choose the mode to be used by the Access Point.

My LAN Connection is : Static IP

STATIC IP ADDRESS LAN CONNECTION TYPE :

Enter the static address information.

IP Address :

Subnet Mask :

Gateway Address :

Helpful Hints..

LAN Settings:
 Also referred as private settings. LAN settings allow you to configure LAN interface of DAP-1160. LAN IP address is private to your internal network and is not visible to Internet. The default IP address 192.168.0.50 with subnet mask as 255.255.255.0.

LAN Connection type:
 The factory default setting is "Static IP" which allows the IP address of the DAP-1160 to be manually configured in accordance to the applied local area network. Enable Dynamic (DHCP) to allow the DHCP host to automatically assign the Access Point an IP address that conforms to the applied local area network.

IP Address:
 The default IP address is

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Dear Client,

Thank you for having chosen our autoclave, we know how to exchange your fidelity, with maximum attention and service definitely corresponding to your expectations.

The autoclave you have chosen is the maximum expression of applied technology. It can be used for dental, medical, and aesthetic applications and in general in all the sectors in which there is the need of performing a safe sterilization of tools and materials of any kind: solid, hollow, and porous loads, both wrapped and non-wrapped.

Before using this autoclave, we invite you to read with maximum attention the user's manual and then keep it in a place accessible to all operators in charge of sterilization.

The autoclave must be used only and exclusively for the uses foreseen by the constructor.

For installation, maintenance and assistance ask exclusively for a technician authorized. We invite you to use and ask for exclusively original spare parts.

CONTENTS

- 1. HOW TO CONSULT THE MANUAL**
 - 1.1 GLOSSARY
 - 1.2 FRONT VIEW / REAR VIEW DESCRIPTION
 - 1.3 DISPLAY DESCRIPTION
- 2. SAFETY**
 - 2.1 SAFETY MARKING
 - 2.2 SAFETY DEVICES
 - 2.3 SAFETY NOTES
 - 2.4 DISPOSAL
- 3. USE AND DESTINATION OF USE OF THE AUTOCLAVE**
- 4. TECHNICAL DATA**
 - 4.1 MECHANIC
 - 4.2 ELECTRIC
 - 4.3 CHAMBER
 - 4.4 CLEAN WATER TANK
 - 4.5 USED WATER TANK
 - 4.6 BACTERIOLOGICAL FILTER
 - 4.7 REMOTE CONTROL CHARACTERISTICS
- 5. ACCESSORIES**
- 6. UNPACKAGING**
- 7. INSTALLATION**
- 8. USE INSTRUCTION**
 - 8.1 TURN ON THE AUTOCLAVE
 - 8.2 AUTOCLAVE SETUP
 - 8.3 CLEAN WATER TANK FILLING
 - 8.4 CHARACTERISTICS OF THE WATER TO BE USED
 - 8.5 LOADING OF THE MATERIALS INTO THE AUTOCLAVE
 - 8.6 LABEL PRINTING FUNCTION
 - 8.7 STARTING THE STERILIZATION CYCLE
 - 8.8 END OF CYCLE
 - 8.9 UNLOADING OF THE STERILIZED MATERIALS
 - 8.10 DISCHARGING USED WATER
 - 8.11 INTERRUPTION OF A STERILIZATION CYCLE
- 9. INTERNET CONNECTION**
- 10. STERILIZATION CYCLES**
 - 10.1 CYCLES DESCRIPTION
 - 10.2 OPERATING CYCLES
 - 10.3 TEST CYCLES

- 10.4 NIGHT / CYCLE
- 10.5 TABLES / CHARTS

11. AUTOCLAVE CHECK TEST

- 11.1 CHEMICAL INTEGRATORS
- 11.2 BIOLOGICAL INTEGRATORS
- 11.3 BOWIE&DICK TEST
- 11.4 HELIX TEST
- 11.5 VACUUM TEST

12. ALARM MESSAGES

13. SOLUTION TO OPERATING PROBLEMS

- 13.1 THE AUTOCLAVE DOES NOT DRY CORRECTLY
- 13.2 THE AUTOCLAVE CHAMBER HAS TURNED WHITE
- 13.3 THE AUTOCLAVE CHAMBER HAS GREEN STAINS
- 13.4 THE STERILIZATION CYCLE
- 13.5 THE AUTOCLAVE DOES NOT RECEIVE CONTROLS
- 13.6 STAINS ON INSTRUMENTS

14. DEMINERALIZER

- 14.1 DEMINERALIZER SETTING
- 14.2 DEMINERALIZER CONNECTION
- 14.3 RESIN AND FILTER CARTRIDGE REPLACEMENT

15. STERILIZATION STANDARD PROCEDURES

16. INSTRUCTIONS FOR THE AUTOCLAVE SHIPMENT

17. SERVICE PROCEDURES

18. WARRANTY

A. MAINTENANCE

B. SOFTWARE VISUALIZATION CYCLES

1**HOW TO CONSULT THE MANUAL****1.1 Glossary**

<i>VACUUM</i>	Pressure below atmospheric pressure
<i>CHAMBER</i>	Cylindrical steel body that, during the cycle, is filled with saturated steam
<i>FILTER</i>	Device for retaining particles
<i>TRAY</i>	Drawer for holding the load to be sterilized
<i>TRAY HOLDER</i>	Structure for holding the trays in the sterilization chamber
<i>DEMINERALIZER</i>	Device for eliminating mineral salts in mains water
<i>HELIX TEST</i>	Test for checking the penetration of steam in a hollow object
<i>BOWIE - DICK TEST</i>	Test for checking the penetration of steam in a porous object
<i>VACUUM TEST</i>	Test for checking that the degree of vacuum in the chamber is maintained
<i>BACTERIOLOGICAL FILTER</i>	Device for retaining particles in suspension in the air larger than 0.3µm.
<i>OVER PRESSURE</i>	Pressure higher than the standard one for the cycle considered
<i>OVER TEMPERATURE</i>	Temperature higher than the standard one for the cycle considered
<i>OSMOSIS FILTER</i>	Device to remove mineral salts from water coming from the water mains, through reverse osmosis
<i>IP ADDRESS</i>	Numerical label (4 digits) identifying univocally a device connected to a computer network (following the IPv4 standard of the IP protocol (Internet Protocol))
<i>SUBNET MASK</i>	Four digits identifying the subnets and the hosts making up the local network (LAN). It is necessary for the correct definition of what IP addresses can be directly accessed and what addresses needs a Gateway.
<i>GATEWAY</i>	Network device whose main aim is to convey network packages outside the local network (LAN), on the demand of the other devices connected to the local network.
<i>DNS</i>	System used for host names to IP addresses resolution and vice versa.

1.2 Front view / rear view description


FIG.1



FIG.2

1	Demineralized water overflow
2	Demineralized water drain tap (<i>demineralizer tap</i>)
3	Used water overflow – condensate drain
4	Used water drain tap
5	Main power supply with fuses
6	Demineralizer connection
7	Ethernet port
8	Spacer
9	Demineralized water pump drain
10	Bacteriological filter
11	RS232 connector
12	<i>B-TEST</i> biological incubator connector
13	Used water drain
14	Memory card slot
B-M1	<i>Multifunction 1</i> button
B-M2	<i>Multifunction 2</i> button
B-M3	<i>Multifunction 3</i> button
B-START	START-STOP button
B-PUMP	Demineralized water filling pump
B-SELECT	Cycle selection button
B-SOS	Assistance Message sending button
B-MAIN	Main switch

1.3 Display description

The display encompasses three informative sections.

1 GRAPHIC INFORMATION SECTION

	Clean water tank empty
	Clean water tank medium level
	Clean water tank full
	Used water tank full
	Water filling pump on
	Connection to the network <i>disabled</i>
	Connection to the network <i>active</i>
	Vocal messages <i>disabled</i>
	Vocal messages <i>active</i>
	Door open
	Door closed

Door locked

Remove control *inactive*

Remove control *active*

Memory card fitted


2 AUTOCLAVE OPERATION INFORMATION SECTION

Information relating to the autoclave status is displayed.


- Cycles selection
- Sterilization cycle phases status
- Display of all of the autoclave working parameters
- Display of all of the autoclave working menus
- Alarms display

3 MULTI-FUNCTION BUTTONS SECTION

The function of the multi-function buttons (BM1 BM2 BM3) depends on the displayed message. Each button corresponds to the function displayed above it. You will find the indication of functions later on in the manual (you will not find the general indication BM1, BM2 or BM3).
For instance:

Indicates that the multi-fuction button corresponding to the displayed arrow has to be pressed .


Unlock	Press to unlock the door
Label	Press to enter the label printing menu
Info	Press to enter the display of the autoclave operating parameters
Setup	Press to enter the autoclave setup menu
Reset	Press the two buttons at the same time to reset the machine after an alarm
Help	Press to enter the help menu
X	Press to cancel a selection
▲	Press to scroll upward the menu or the alphanumerical values
▼	Press to scroll downward the menu or the alphanumerical values
—	Press to move to the next alphanumerical value
✓	Press to confirm a selection

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DISCONNECT THE POWER BEFORE
REMOVING THE LID

WARNING
DISCONNECT THE
MAINS SUPPLY BEFORE
REMOVING THIS COVER

ATTENZIONE
ATTENTION
ATTENTION
ACHTUNG

HIGH TEMPERATURE

ALTA TEMPERATURA
HIGH TEMPERATURES
TEMPERATURES ELEVÉES
HOHE TEMPERATUR

EARTH CONNECTION

2.2

Safety Devices

The following safety devices are installed:

- Safety valve set at 2.4 bar 0/+10%
- Electromagnetic lock to prevent the door from opening while the cycle is running
- Resistance over temperature thermostats
- Safety pressure switch working on the automatic door lock

2.3

Safety Notes

- The manufacturer is liable for the marketed product in accordance with current regulations. The manufacturer's liability will expire when operations are carried out on the device, or a part of it, by unskilled personnel or using non-original spare parts.
- There should be no potential risk of explosion and/or fire in the room where the autoclave is installed.
- The autoclave should be installed in a special well-ventilated room.
- The incubator's connector (FIG.2-pos.12), should be used only for the B-TEST.

Since it is a 230 V a.c. connection independent from the main switch, do not touch the connector with wet or moist hands and protect it against possible water splashes.

2.4 Disposal

This product is subject to Directive 2002/96/EC of the European Parliament and the Council of the European Union on Waste of Electrical and Electronic Equipment (WEEE) and, in jurisdictions adopting that Directive, is marked as being put on the market after August 13, 2005, and should not be disposed of as unsorted public waste. Please utilise your local WEEE collection facilities in the disposition and otherwise observe all applicable requirements.

3 USE AND DESTINATION OF USE OF THE AUTOCLAVE

The autoclave is able to sterilize the three types of load provided for by the standard EN13060, especially:

METAL OR SOLID MATERIALS	Instruments with no cavities and no obstacles to the penetration of steam
POROUS OBJECTS	Simple or composite materials that can absorb fluids (fabrics, gowns, surgical gauzes, dressings, etc)
HOLLOW OBJECTS	Materials or devices with cavities, obstructions, etc. These are subdivided into two types, classified according to the length and diameter of the cavity. Approximately: TYPE B: cannulas, tubes or devices with large passages. TYPE A: turbines, hand pieces and devices with blind or small holes

This device has been certified for the sterilization of the following materials *:

	AURA B	AURA B 24
Wrapped and unwrapped solid materials for a maximum load of	max kg. 5	max kg. 7
Type A/B hollow materials for a maximum load of	max kg. 5	max kg. 7
Porous materials for a maximum load of	max kg. 1,5	max kg. 2

* Only for European Countries

The autoclave shall only have to be used for the sterilization of tools and materials being compatible with the steam sterilization system. Always make sure that the loads that need to undergo sterilization can stand the temperatures scheduled for the selected cycle.

4
TECHNICAL DATA
4.1 Mechanical Data

	AURA B	AURA B 24
Working temperature	+5°C ÷ +40°C	
Maximum altitude	2.000 m (s.l.m.)	
MAX relative humidity at 30°C	80%	
MAX relative humidity at 40°C	50%	
Dimensions of space occupied (L x H x P)	504 X 400 X 640 mm	504 X 400 X 730 mm
Space occupied with open door	300 mm	
Weight (tank full)	70 kg	74 kg
Weight of area of support	2058 N/m ²	2058 N/m ²
Potential sound level	< 70 db A	

4.2 Electrical Data

	AURA B	AURA B 24
Power voltage	230 V a.c. +/-10 % single phase	
Power	1,7 kW	2.2 kW
Frequency	50 / 60 Hz	
Power cord	2 + 1 x 1mm ²	
Fuses	5 x 20 - 10A	
Heat transmitted	3.6 E ⁶ J / ora	

4.3 Chamber

	AURA B	AURA 24 B
MAX working pressure	2.4 bar (relativi)	
MAX empty	- 0.9 bar (relativi)	
MAX Temperature	138 °C	
Material	Inox AISI 304	
Size	Ø 245 x 318 mm	Ø 245 x 500 mm

4.4 Clean Water tank

	AURA B	AURA B 24	
		With Demineralizer	Without Demineralizer
Volume	4,5 l		
Usable cycles	4	Read demineralizer instruction	2
Material	polyethylene		

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





Tray extraction and door adjustment wrench

Use for extract the trays (FIG.11) and for door adjustment.

Picture

FIG.10


FIG.11

Envelope standard

1

Chamber and door gasket clearing sponge

Use to clean sterilization chamber and door gasket (par.A)


FIG.12

Envelope standard

1

Memory card

To be used to memorize the autoclave cycles (for the best compatibility with the device, it is advisable to always use the original memory card).

WARNING: The card contains the software reading log cycle - Carry out the rescue and installed on the PC before the commissioning of the autoclave (the instruction manual is on the card – par.B)

FIG.13

Envelope standard

1

Water filling pipe with filter

Use to charge water with pump (*par.8.3*).

FIG.14

Envelope standard

1

Water discharge pipe

Use to discharge used water (*par.8.9*). Connect the metal fitting into the demineralized water charge port (*FIG.2–pos.13*)


FIG.15

Envelope standard

1

Water discharge pipe / condensate drain

Connect one pipe end into the used water overflow - condensate drain pipe fitting (*FIG.1–pos.3*), and put the other end in a tank.


FIG.16

Envelope standard

1

Rear spacer

Put the spacer in the autoclave's back panel (*FIG.1–pos.8*) It's necessary for guarantee a good ventilation if you place the autoclave near a wall.


FIG.17

Envelope standard

1

Pipes for discharge utilities

1- Pipe for clean water overflow: Connect one pipe end into the water overflow pipe fitting (*FIG. 1–pos. 1*), and put the other end in a tank or in the discharge (demineralizer version).

2- Pipe for discharging used water: Connect one pipe end at the pipe fitting (*FIG.20*) and screw it into the used water tap (*FIG. 1–pos.4*), then put the other end in a tank or in the discharge.


FIG.18

Envelope standard

2

Pipe fitting for back discharge used water

Connect the pipe for discharging used water into the pipe fitting for back discharge used water (*FIG. 19*) and screw into the used water tap (*FIG. 1–pos.4*).


FIG.19

Envelope standard

1

Power supply cord

Take the power supply cord provided and insert the female plug (*FIG.21–pos.B*) in the socket of the back panel of the autoclave (*FIG.1– pos.5*) Then insert the male plug (*FIG.21–pos.A*) in the electric plug of the system.



FIG.20

Envelope standard

1

Remote control / Remote Base

The autoclave is fitted with remote control (remote base) to select and start cycles and to have remote information on the end of the cycle.

6

UNPACKING

The autoclave is shipped in a suitable package to be transported and moved easily and to protect its contents.

The package must not be jostled, it must be moved with care avoiding letting it turn over or fall from a height superior to 16 cm.

In case autonomous handling means are not available handle the packaged autoclave always in two persons.

Autoclave is placed on wooden pallet and packed with corrugated cardboard internal and external application.

For remove the autoclave to its pack please remove first the corrugated cardboard.

For lift the autoclave please use the belts.

The autoclave shall have to be handled by at least two people and by using belts only.

Do not lift the autoclave from the inferior part of the door or command panel, this incorrect operation can create problems of a mechanical nature.

ATTENTION: Always conserve original packaging.

Inside the package you will find:

- *USER'S MANUAL*: to be read with attention and kept in a place available to all operators assigned to sterilization.
- *CERTIFICATION*: which must be conserved.
- *GUARANTEE LABEL*: which must be completed and send.


FIG.21

	AURA B	AURA B 24
A	640 mm	730 mm
B	50 mm	
C	400 mm	
D	504 mm	
E	270 mm	


- 1 – Install the autoclave in environments suitable to sterilization.
- 2 – The premises have to be suitably lit and ventilated, in compliance with regulations in force.
- 3 – Install the autoclave away from heat sources and water splashes.
- 4 – Position the autoclave on a plane suitable to bear its weight (min.80 kg) and of suitable dimensions.
- 5 – The surface supporting the autoclave must be perfectly horizontal.
- 6 – Place the autoclave at a height allowing easy intervention by the operator for inspection and clearing of the whole sterilization chamber.

- 7 – Open the door of the autoclave and remove all packages containing the single accessories from inside the sterilization chamber.
- 8 – Leave inside the sterilization chamber only the tray carrier with the trays. All other accessories should be positioned in a separate space available to operators.
- 9 – Never place on the autoclave newspapers, trays, liquid containers, etc...
- 10 – Do not lean on the door when it is open.
- 11 – Leave a space of at least 5 cm in the rear of the autoclave, using the spacer (*FIG. 1–pos.8 / FIG. 17*) and in the sides of the unit to ensure the ventilation required.
- 12 – Connect the overflow discharge hoses in the rear panel (*FIG. 1–pos. 1-3*). Water can be discharged either in a tank or directly into the tub discharge.
- 13 – Always make sure the electrical system to which the autoclave is to be connected is in conformity with the legislation in force and sized to suit the specifications of the said device.
- 14 – Connect the power supply cord to the socket on the rear panel of the autoclave (*FIG. 1–pos.5*)
- 15 – Connect the power plug to the system and make sure the power supply is 230 Vac-2000Va

IMPORTANT:

Do not use extensions, reducers or adapters for connection as this could cause micro interruptions with consequent generation of alarm signals.

- 16 – In the end, to prepare the tool for the network connection, switch it off, fit the network cable RJ45 (*not provided with the tool*) into the port in the rear panel (*FIG. 2–pos.7*). Connect the other connector to your network connection device (*Hub, Switch, Router or PC*)

8.1 Turn on the autoclave (after the first switch on)

Press the main switch **B-MAIN** (FIG.2)

After the display of the logo, the autoclave checks the memory and the connections. The result of the checks is shown on the display (FIG.22).

FIG.22

Once checks are over the autoclave shifts to the operating setup (FIG.23).


FIG.23

Open the door and wait for some seconds until an acoustic signal informs that the automatic barometric alignment values have been acquired; at the same time the display shows the message DOOR OPEN. The autoclave is ready for use.

IMPORTANT:

Never select an item before the acoustic signal, the autoclave will not accept the preset selection.

8.2 Autoclave setup

From the first page press **Setup** to enter the autoclave setup menu.

FIG.24

LANGUAGE:ENGLISH

Press the central button to change the language of menus and vocal messages

Press the arrow to shift to the next item

DATE AND TIME

Press the central button to enter the date/time setup.

00:00

00/00/1900

When the box flashes press the arrows to select the desired value.

00:00

00/00/1900

Press once more the central button to move among boxes and the arrows to select a value.

Go on until the last value is selected.

Press the central button for the last time to display the final selection page.

Press to confirm selected values and go back to the setup menu

Press to restart the procedure

Press to cancel selected values and go back to the setup menu

Press the arrow to shift to the next item

LAN SETUP

Press the central button to enter the LAN network setup menu

- **INDIRIZZO IP** : Assign the four digits following your network configuration. Carry out the same procedure as the one of the date/time menu to select and validate the entered data and confirm the entered value. The default value is 192.168.0.1
- **SUBNET MASK** : Assign the four digits following your network configuration. Carry out the same procedure as the one of the date/time menu to select and validate the entered data and confirm the entered value. In most cases the default value 255.255.255.0 is compatible with your configuration.
- **GATEWAY** : Assign the four digits following your network configuration. Carry out the same procedure as the one of the date/time menu to select and validate the entered data and confirm the entered value. The default value is 192.168.0.254
- **SERVER DNS** : Assign the four digits following your network configuration. Carry out the same procedure as the one of the date/time menu to select and validate the entered data and confirm the entered value. The default value 208.67.222.222 corresponds to an existing global DNS server. We suggest to keep the preset value.
- **VALORI DI DEFAULT**: Resets the preset values of the previous parameters of LAN SETUP menu.

CONSUMABLES AND SPARE PARTS

Press the central button to enter the menu for consumables and spare parts orders.

REMOTE CONTROL INSTALLATION

- Place the batteries (do not reverse correct polarity) and switch on the remote control through the switch.
- To connect the remote control to the autoclave enter the menu. Hold down for 5 seconds one of the two buttons of the remote control (a sound and a vibration will confirm that the operation has been successful).
- Now the remote control is active.

INFORMATION MESSAGES

Press the central button to enable/disable the information messages

OSMOSIS CONTROL

The menu allows the activation of the feeding system by the osmosis device (optional). When the system is activated the charging pump is disabled.

The cycles carried out since the last filters replacement are also shown.

The counter cannot be reset when filters are replaced.

See chapter 14 for further details.

VOCAL MESSAGES

Press the central button to enable/disable vocal messages

SERVICE MODE

Access to the service menu (after password entry). This mode is reserved exclusively for the setup by an authorised engineer.

The manufacturer shall not be held responsible for any tampering with or injury to unauthorised staff.

8.3 Clean water tank filling

Connect the water filling hose with filter (*FIG.14*), and fit it into the front coupler of the autoclave (*FIG.2- pos.9*).

Fit the other end of the hose with filter inside the tank of demineralized or distilled water.


Now press the button **B-PUMP** (*FIG.2*) to start the water filling pump. If the maximum level is not achieved within 180 seconds the pump stops automatically and it will be necessary to press button **B-PUMP** again to complete the tank filling. The pump stops automatically when the maximum level is achieved.

If the autoclave is mounted with the demineralization system comply with the instructions shown in the package of this device and those of chapter 14.

8.4 Characteristics of the water to be used

TABLE SHOWING THE QUALITY LEVELS LAID BY THE UNI EN 13060: 2009

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8.5 Loading of the materials in the autoclave**FIG.25**




Place the tray carrier into the sterilization chamber in horizontal position by arranging the materials to be sterilized on the provided trays, as follows:

- Do not superpose the materials
- Arrange the wrapped materials with the paper side facing upwards
- Never bring the materials into contact with the sterilization chamber or the autoclave door
- Put scissors and dental forceps with blades open

After loading the materials, close the door of the autoclave. The display will show the icon and the message DOOR CLOSED.

8.6 Label printing function

After the printer installation (*optional accessory*), from the initial page by pressing the button **Label** the menu for the printing of the traceability adhesive labels to be stuck on the wrap before starting sterilization is entered.

	<p>When the box flashes press the arrows to select the number of labels to be printed.</p>
	<p>Press the central button to move to the next selection</p>
	<p>When the box flashes press the arrows to select the number of months before expiry. Press the central button for the last time to display the final selection page.</p>

8.7 Starting the sterilization cycle

Select the sterilization cycle most suitable to the loaded materials by pressing the button **B-SELECT** (FIG.2) Each single pressure changes the program selection to the next one. Programs are displayed (FIG.23).

After the program selection, start the cycle by pressing the button **B-START**(FIG.2). The door is locked automatically and the cycle starts.

During the cycle the display shows all parameters and information related to the cycle in progress. In this setup the display shows: the type of cycle, the cycle status, the time remaining before the end of cycle (for the Vacuum Test it identifies the whole cycle, while for all the other cycles it identifies the sterilization phase plus that of drying), the number of cycles done by the machine and the button **Info** allowing access to the list of working parameters (FIG.26).


FIG.26

8.8 End of cycle

An acoustic signal and a displayed image warn that the sterilization cycle has been completed. Press any button to continue and unlock the door by pressing one of the buttons **Unlock**. If there is pressure inside the chamber the button will not activate the release device. When the door is unlocked pull the door handle and open.

8.9 Unloading the sterilized materials

Wear personal protective equipment in compliance with regulations on safety and hygiene at work. Extract the trays by using the spanner provided (FIG.10), let rest the instruments and store them in environments which are not exposed to contamination.

8.10 Discharging used water

When the icon of used water turns on the tank has to be discharged, otherwise the autoclave operation is inhibited.

Fit the hose (*FIG.16*) into the water discharge plug on the front panel of the autoclave (*FIG.2-pos.13*), at the same time put the other end of the hose into a container; the water will drop into such container.

IMPORTANT:

A - The hose located in the receiving container, must never touch or be immersed in the discharged water, otherwise sucking up will occur.

B – Always wait until the discharge water has been completely drained. The used water maximum level LED turns off when some water still remains into the tank, consequently do not consider it as a reference for this operation.

At the end of the discharge press the ring nut on the plug and remove the hose.

8.11 Interruption of a sterilization cycle

A sterilization cycle can be voluntarily interrupted by pressing the button **B-START** for at least 2 seconds.

The autoclave will send out a beep, align with the atmospheric pressure and the display will show the alarm message A001 (*Interrupted cycle*).

To reset the alarm hold down at the same time the buttons below the bar indicated with **Reset** (*FIG.27*).

If an alarm message is displayed (suffix "A") the cycle has to be considered as **FAILED**: all loading and sterilization operations have to be repeated.


FIG.27

9**INTERNET CONNECTION**

The autoclave AURA B requires internet connection to be able to use the services created for it.

To connect the autoclave to your network, TECNO-GAZ S.p.A proposes the following connecting systems to suit your needs:

- *Direct connection to your LAN with network cable (recommended);*
- *Connecting to your LAN through WI-FI bridge (if the network connection with cable is not feasible);*
- *Replacing your ADSL modem with a router ADSL WIFI and WIFI bridge to link autoclave (if you don't have a LAN);*

9.1 Direct connection to your LAN with network cable RJ-45


Connect the network cable RJ-45 autoclave (*FIG.1-pos.7*) and the other end to your router (or switch or hub); To configure your autoclave you need the following information (already in your possession or alternatively contact your system administrator):

- Free IP address belongs to your LAN;
- Gateway (usually the IP address of your router);

Enter in the user menu of your autoclave, revenues in the **LAN SETUP** (see par.8.2) and insert the above parameters in their respective fields. After the procedure, turn off and switch on again the autoclave.

If the parameters you entered are correct, during the power phase (*FIG.22*), the autoclave should verify that the server is operational TECNO-GAZ (voice "Online Server"). Otherwise check the data entered in the preceding paragraph. If the problem persists, contact your system administrator. Also, if the codes have been entered correctly (*Customer ID*) and reseller (*Reseller ID*), the autoclave will show a confirmation screen of the activated Internet services associated with your autoclave.

9.2 Connecting to the LAN by a WI-FI bridge

You can create a WI-FI bridge between autoclave and nearest outlet. To proceed to the configuration of this system is necessary to access the LAN configuration of the three devices:

- **Configuring Access-point WI-FI**
- Connect the 'Access Point-WI-FI LAN via network cable to your computer;
- Connect the power cable to the device;
- Wait 30 seconds that LED on the device light up;
- Open your Internet browser and type in the address: 192.168.0.254;

LOGIN

Log in to the Access Point

User Name :

Password :

Log In

- Insert the following parameters (being careful Case/Lowercase):
 - Username: admin
 - Password: AuraBPwd
- Press “Log In”;
- Go to the section: ”SETUP -> Lan Setup”;

Product Page: DAP-1160
 Firmware Version: v1.20

DAP-1160	SETUP	ADVANCED	MAINTENANCE	STATUS	HELP
----------	-------	----------	-------------	--------	------

Setup Wizard
 Wireless Setup
 LAN Setup
 Logout
 Reboot

NETWORK SETTINGS :
 Use this section to configure the internal network settings of your AP and also to configure the built-in DHCP Server to assign IP addresses to the computers on your network. The IP Address that is configured here is the IP Address that you use to access the Web-based management interface. If you change the IP Address here, you may need to adjust your PC's network settings to access the network again.
 Save Settings Don't Save Settings

LAN CONNECTION TYPE :
 Choose the mode to be used by the Access Point.
 My LAN Connection is : Static IP

STATIC IP ADDRESS LAN CONNECTION TYPE :
 Enter the static address information.
 IP Address : 192.168.2.253
 Subnet Mask : 255.255.255.0
 Gateway Address : 192.168.2.1

Helpful Hints..
LAN Settings:
 Also referred as private settings. LAN settings allow you to configure LAN interface of DAP-1160. LAN IP address is private to your internal network and is not visible to Internet. The default IP address 192.168.0.50 with subnet mask as 255.255.255.0.
LAN Connection type:
 The factory default setting is "Static IP" which allows the IP address of the DAP-1160 to be manually configured in accordance to the applied local area network. Enable Dynamic (DHCP) to allow the DHCP host to automatically assign the Access Point an IP address that conforms to the applied local area network.
IP Address:
 The default IP address is

- Enter the followings parameters:
 - IP Address;
 - Gateway Address;
- Press “Save Settings”;
- Close your browser;
- Disconnect the Access-Point WI-FI from your computer;

- **WI-FI client configuration:**

Proceed as described in the previous case to assign an IP address and a gateway to the WI-FI client, Disconnect WI-FI client from your computer and proceed to configure the system;


- **WI-FI bridge connect:**

- Connect the Access-Point WI-FI to the available outlet;
- Connect the power cable to the WI-FI Access-Point;
- Wait 30 seconds that the WI-FI LED will activate;
- Connect the WI-FI Client to the autoclave with RJ-45 cable;
- Connect power cable to the WI-FI Client;
- Wait 30 seconds that WI-FI LED on the device light up;

At this point the WI-FI bridge is already established between the two devices, turn on your autoclave and proceed to configure the same as in the previous section (*par. 9.1*).

9.3

Replacing your ADSL modem with a router ADSL WIFI and WIFI bridge to link autoclave

If your internet system connection using an ADSL modem (USB or Ethernet) providing access to one computer, you can replace your modem with a router WI-FI ADSL. In this way you will have the opportunity to share the internet connection to your PC and to your autoclave.

- **Configuring WI-FI Router ADSL:**

- Unplug your ADSL modem and connect the phone line to your new router by inserting the cable into the port labeled "DSL".
- Connect the RJ-45 network cable to "LAN 1" port and to your PC.
- Connect the power cable and turn on your device.
- After about 30 seconds open your favorite Internet browser and type in the address: 192.168.0.250
- Insert the following parameters (being careful Case/Lowercase):
 - Username: admin
 - Password: AuraBPwd

LOGIN

Log in to the Access Point

User Name :

Password :

- Press "Log In"
- Go to the section: "SETUP -> ADSL Setup" and enable the flag "Manual Setup"

DSL-2640B // **SETUP** **ADVANCED** **MAINTENANCE** **STATUS**

ADSL Setup
Wireless Setup
LAN Setup
Time and Date
Parental Control
Logout

ADSL SETUP

If you are configuring this device for the first time, D-Link recommends that you click the Setup Wizard button, and follow the instructions on screen. If you wish to modify or configure the ADSL settings manually, tick Manual Setup to enable the ADSL Connection Setup.

☒ **Manual Setup**

- Select the type of ADSL connection you have available (in most cases PPPoE / PPPoA)
- Introduce the following parameters (provided by your ISP) for your ADSL connection:
 - Username;
 - Password;
 - Service Name (es MyADSL);

- Press “Save Settings”;

MANUAL ADSL CONNECTION SETUP

Please select the appropriate option to connect to your ISP.

☒ **PPPoE/PPPoA**
Choose this option if your ISP uses PPPoE/PPPoA.(For most DSL users)

☐ **Dynamic IP Address**
Choose this option if your ISP uses Dynamic IP Address over DSL.

☐ **Static IP Address**
Choose this option if your ISP uses Static IP assignments.

☐ **Bridging**
Choose this option if your ISP uses Bridging.

PPPOE/PPPOA INTERNET CONNECTION TYPE :

Enter the information provided by your Internet Service Provider (ISP).

Username :

Password :

Service Name :

Connection Type :

PPPoE LLC

MTU :

1492

Idle Time Out :

0

Minutes (0 = Always On)

VPI :

8

VCI :

35

Save Settings

Connect

Cancel

- Go to the section “SETUP -> Lan Setup”
 - Introduces (possibly) a new IP address per device;
 - Press “Save Settings”;

Product: DSL-2640B

Firmware Version: EU_4.00 Hardware Version: B2

DSL-2640B

SETUP

ADVANCED

MAINTENANCE

STATUS

HELP

ADSL Setup
 Wireless Setup
LAN Setup
 Time and Date
 Parental Control
 Logout

Internet
 Offline

English

Reboot

LAN SETUP

This section allows you to configure the local network settings of your router. Please note that this section is optional and you should not need to change any of the settings here to get your network up and running.

ROUTER SETTINGS

Use this section to configure the local network settings of your router. The IP Address that is configured here is the IP Address that you use to access the Web-based management interface. If you change the IP Address here, you may need to adjust your PC's network settings to access the network again.

Router IP Address :

192.168.1.250

Subnet Mask :

255.255.255.0

DHCP SERVER SETTINGS (OPTIONAL)

Use this section to configure the built-in DHCP Server to assign IP addresses to the computers on your network.

Enable DHCP Server :

☐

DHCP IP Address Range :

to

DHCP Lease Time :

(hours)

Helpful Hints...

If you already have a DHCP server on your network or are using static IP addresses on all the devices on your network, uncheck Enable DHCP Server to disable this feature.

If you have devices on your network that should always have fixed IP addresses, add a Static DHCP for each such device.

More...

- **WI-FI Client configuration:**

Proceed as described in the previous case to assign an IP address and Gateway (on your new ADSL router) to the client WI-FI; disconnect WI-FI Client from your computer and proceed to configure the system;

- **Connecting WI-FI bridge:**

- Connect the WI-FI Router to the ADSL;
- Connect the power cable to the WI-FI router;
- Switch on the device and wait 30 seconds that the WI-FI LED are activated;
- Connect the WI-FI Client to the your autoclave with RJ-45 cable;
- Connect the power supply to WI-FI Client;
- Wait 30 seconds that the WI-FI LED will activate;

At this point the WI-FI bridge is already established between the two devices, turn on your autoclave and proceed to configure the same as in the previous section.

10**STERILIZATION CYCLES****10.1 Descriptions cycles**

The autoclave has three series of cycles:

- Operation cycles
- Night cycles
- Test cycles

10.2 Operation cycles

All operation cycles have a system of fractioned vacuum, which therefore permits sterilization of materials that are hollow, porous, solid, both free and packaged.

Temperatures can be selected from 121°C – 134°C.

Normally the cycles of 121°C are used for thermoplastics or sensitive materials, while the 134°C cycles are used for all other materials.

In all cases always follow the indications given by the manufacturer of the instruments or devices to be sterilized.

There are also two special cycles:

- Prion cycle, studied for the Crutzfeldt-Jacobs disease (mad cow disease)
- Flash cycle created for rapid sterilization of unpackaged instruments and devices.

10.3 Cicli test

The available test cycles are the following:

- Helix e B&D
- Vacuum test

10.4 Night cycle

The autoclave is provided with a special economizer device.

This device allows the execution of sterilization cycles without the presence of an operator.

At the cycle end, if the door is not opened, the autoclave stabilizes until the intervention of the operator.


10.5 Tables / Charts
STERILIZATION TABLE

CYCLES	EXPOSED TIME T4 (Minuts)	DRYING TIME T5 (Minuts)		PRESSURE WORK MIN AND MAX (relative bar)	TEMPERATURE WORK MIN AND MAX (°C)	
		18L	24L			
121°C (wrapped and unwrapped)	18	12	13	1.04 ÷ 1.30	121 ÷ 125	OPERATIVE CYCLE
134°C (wrapped and unwrapped)	5	12	13	2.02 ÷ 2.38	134 ÷ 138	
134°C Flash	4	3	4	2.02 ÷ 2.38	134 ÷ 138	
134°C Prion	20	12	13	2.02 ÷ 2.38	134 ÷ 138	
134°C Helix / Bowie&Dick	3.5	4		2.02 ÷ 2.38	134 ÷ 138	TEST CYCLE
Vacuum				-0.9		

The pre-heating time and fractionation of the vacuum can vary depending on the conditions of the autoclave from 25 - 35 minutes that will add up with the times of the cycles shown in table. All sterilization cycles have 3 stages of vacuum.

TABLE OF TESTS TYPE AS RULE EN13060:2009	
Type test	Operative cycle
Dynamic sterilizer chambre pressure	X
Air Leakage	X
Empty chamber	X
Solid load	X
Small porous items	X
Small porous loads	X
Full porous load	X
Hollow load B	X
Hollow load A	X
Multiple wrapping	X
Dryness, solid load	X
Dryness, porous load	X

* in compliance with all applicable clauses of standard

CYCLE DIAGRAM

T1, T2, T3
pre heater fase and
vacuum fases

T4
sterilization

T5
drying

FIG.28

11

AUTOClave VALIDATION TESTS

11.1 Chemical integrator


These tests are based on the properties of indicator dyes that are able to modify their own aspect if they are adequately exposed to heat and pressure according to given times, since they interact with all the cycle parameters (pressure, temperature, time). If the test is successfully completed, the correct stability is validated.

The chemical integrators (*spare part 200/S e 215-S*) are to be put into the sterilization chamber before starting the cycle, even with the load already in.

The integrators can either be unwrapped and positioned into the trays or they can be wrapped into the sterilization bags together with the materials to be sterilised.

In any case it is advisable to wrap the integrators where there are many operators. As a matter of fact by doing so the sterilization of every single tool will be certified.

*If the colour change remains within the white area, the test **failed***

*if the colour change enters the **SAFE** window, the test **succeeded***

There is no preset testing interval, apart from areas to which specific regional regulations apply. In any case all the chemical integrators are the cheapest and quickest available tests, hence in order to constantly monitor the autoclave performance, **it is advisable to use them on all of the cycles or at least once a day.**

11.2 Biological indicator

The test (*spare part 262-S*) is used to prove the autoclave capability as to the complete destruction of all microorganisms. It is represented by a standardized spore preparation, having biological characteristics, high resistance to heat, which is to be deemed an absolutely safe control tool (ATCC 7953). It consists of spore ampoules. It is non pathogenic, non toxic and non pyrogenic.

The ampoules shall have to be put into the sterilization chamber even when loaded.

Carry out the cycle and then take out the ampoule by taking care since it is hot and pressurized.

Leave it cool down for about 10 minutes and then activate it by using the supplied crusher by taking care so as to keep the ampoule always upright. The chemical indicator placed on the label of the ampoule will have changed its colour from blue into black.

Subsequently incubate the ampoule into a biological incubator at 57°C for 48 hours, together with an unprocessed test ampoule, that has anyway be activated through the relevant crusher.

After 48 hours take the processed ampoule out of the incubator and assess the test result. If the ampoule has changed its colour to yellow this means that the autoclave has failed the test and there is bacterial growth. If the ampoule remains purple, this means that there is no bacterial growth and hence the test succeeded. The test ampoule of course will always change to yellow since it has not been processed and will only be used as term of comparison.

At the end of the test dispose of the ampoule in the solid urban waste. It is anyway advisable to perform another sterilization cycle at 121°C on the ampoule.

There is no preset testing interval, apart from areas to which specific regional regulations apply. In any case the biological tests are the more concrete ones. Hence in order to constantly monitor the autoclave performance, **it is advisable to perform the test at least once every 90 days.**

11.3 | Bowie-Dick test

It is a physical test (*spare part 268/S*) enabling to verify steam penetration capability into porous loads. For the test a “test pack” is used which is standardized and complying with the technical standards being currently in force.

Perform the test with empty chamber. The Bowie-Dick test pack shall have to be positioned on the autoclave central tray.

Perform the cycle, indicated in the display of the autoclave, and at the end remove the pack, open it and check the sheet colour change through the chemical indicator placed inside.

The assessment is simple and rapid. If the colour change turns out to be uniform (as shown in the figure) the test perfectly succeeded. On the contrary the test failed, hence the autoclave cannot properly sterilize porous bodies.

There is no preset testing interval, apart from areas to which specific regional regulations apply. Hence in order to constantly monitor the autoclave performance, **it is advisable to perform the test at least once every 30 days.**


11.4 Helix test

It is a physical test (*spare part 267-S*) enabling to verify steam penetration capability into hollow loads. For the test a standardized system is used which complies with the technical standards in force.


Put the test device into the empty chamber.

Fit the suitable indicator strip into the test device end cap, then position the test device on the central tray of the autoclave.

Perform the suitable cycle, indicated in the autoclave display, then remove the test device, open the cap and check the indicator strip colour change.

PERIODICITY There is no preset testing interval, apart from areas to which specific regional regulations apply. Hence in order to constantly monitor the autoclave performance, **it is advisable to perform the test at least once every 30 days.**

11.5 Vacuum test

It is a chamber leak test or a vacuum leak test. This test aims at verifying that during the cycle no air infiltration takes place through the chamber seals (gaskets, valves, etc...)

USE Perform the test with empty chamber.

Select the suitable cycle indicated on the autoclave display. The cycle is enabled.

The autoclave automatically performs the cycle by complying with precise technical procedures. Then the final result is printed out through the printer connection.

RESPONSE The test result is immediate and is printed by the printer which prints out all the values relating to the performed cycle as well as the final assessment.

PERIODICITY There is no preset testing interval, apart from areas to which specific regional regulations apply. Hence in order to constantly monitor the autoclave performance, **it is advisable to perform the test at least once per week.**

**WE RECOMMEND TO USE ONLY TEST PROPOSED
BY TECNO - GAZ S.p.A**

12

ALARM MESSAGES

Alarm messages are highlighted by an alphanumerical code made of a letter and 3 digits.

The prefix "A" refers to alarms, anomalies of the autoclave; if after the remedy has been tried the problem remains phone assistance has to be required.

If an alarm message is displayed (suffix "A") the cycle has to be considered as **FAILED**: all loading and sterilization operations have to be repeated.

To reset alarms and errors hold down at the same time the buttons under the bar indicated with **Reset**.

ERROR	CAUSE
A 001	Cycle interrupted by the user
A 101	Vacuum not achieved in 10 min.
A 111	Vacuum not maintained on the first phase of VACUUM TEST
A 121	Vacuum not maintained on the second phase of VACUUM TEST
A131	During the warm-up phases the tool has not charged the correct amount of water
A132	Error in the flow-meter operation
A133	Pressure above the allowable limit during the water recall for the warm-up phases
A 200	Error in the control of the SV operation
A 400	Error in the locking operation
A 401	Error in the hatch operation
A 402	Error in the hatch operation
A 403	Error in the hatch operation
A 404	Error in the hatch operation
A 405	Error in the hatch operation
A 551	Pressure out of limit
A 637	Error in the access to the SD card
A 651	Reading of the T1 probe during sterilization above the maximum limit
A 653	Reading of the T2 probe during sterilization above the maximum limit
A 661	Error in the probe reading
A 662	Error in the probe reading
A 701	Error for failed pressure achievement during the first two phases of the warm-up
A 711	Error for failed pressure achievement during the third phase of the warm-up
A 751	Reading of the T1 probe during sterilization below the minimum limit

A 753	Reading of the T2 probe during sterilization below the minimum limit
A 781	Temperature of the 121°C cycle out of the maximum limit
A 782	Temperature of the 134°C cycle out of the maximum limit
A 801	Errors of out of maximum time during the first phases of discharge
A 811	Errors of out of maximum time during the last phase of discharge
A 901	Cycle interrupted due to power supply failure

For all alarms it is advisable to reset the machine and restart the cycle.

If the problem remains ask for the phone service at MEDILINE ITALIA srl +39 0522 - 94 29 97 or through the button **SOS** (FIG.2). In this way the device, if connected to the network, will send a service request to our center. You will be contacted by one of our engineers as soon as possible.

13 SOLUTION TO OPERATING PROBLEMS

In most cases alarms or errors are caused by lack of attention or lack of familiarity with some technical and operational aspects.

You will find below the list of some anomalies with associated solutions.

13.1 The autoclave does not dry materials correctly

- Replace the bacteriological filter with a new original one.
- Non original trays, without holes or with different holes made of different material have been used. It is advisable to use original trays only.
- Instruments have not been arranged correctly. Carefully comply with instructions of par.8.4

13.2 The autoclave chamber has turned white

- Change immediately the type of water used, use exclusively demineralized or distilled water, as indicated in the previous chapters and then clean the chamber.
- The whitish colour may be a consequence of the evaporation of organic materials from the instruments. Instruments should undergo a more suitable and deeper cleaning .
- Check the demineralizer installed.

13.3 The autoclave chamber has bluish-green stains

- Instruments have not been cleaned correctly after the cleaning, rinse more carefully, if the stains remain apply to the phone technical service .

13.4 The sterilization cycle interrupts without any patent reason

- Check whether the autoclave is connected to the power mains through extensions, reducers or adapters, if so remove such accessories and connect the autoclave directly to the power socket.

13.5 The autoclave does not receive controls

- The autoclave is carrying out the automatic barometric alignment, wait for the double acoustic signal after the door opening , then set functions.
- The demineralized water tank is empty, the Led indicating the minimum level is ON, fill the tank with pure water.
- The used water tank is full, the Led indicating the maximum level is ON, discharge used water .

13.6 Stains on instruments

- Instruments get yellow due to the deposit of residual chemical fluid which has fixed on the instruments due to heating. The rinsing was not suitable.
- The sterilization chamber has yellow stains. Some instruments with traces of chemical fluid have been loaded into the chamber, such fluid has fixed due to heating. The rinsing was not adequate.
- Instruments have whitish stains, the rinsing has been done with very limy water and the instruments have not been carefully dried. For the final rinsing we suggest the use of demineralized water and a careful drying of the instruments
- Instruments have blackened since the material they are made of contain a high percentage of carbon.

14**DEMINERALIZER****14.1 Demineralizer setting**

The autoclave can also charge demineralized water through the *PURE OSMOSIS SYSTEM PLUS* (Optional accessory).

This application is used to remove mineral salts from the water coming from the water network before filling the autoclave autoclave.

14.2 Connection of the demineralizer

You will find below details on the specific connection to the autoclave.

Turn off the autoclave, if ON (*FIG.2*)

- Close the cock upstream the demineralizer
- Install the demineralizer as indicated in its manual;
- Wrap the male thread of the hose coupler with teflon or another component ensuring water tightness;
- Screw the hose coupler on the female thread of the clean water discharge (*FIG.1–pos.2*)
- Insert the hose from the demineralizer into the hose coupler and screw it to the autoclave;
- Insert the feeding coupler of the demineralizer into the port (*FIG.1–pos.6*) on the rear panel of the autoclave;
- Open the valve upstream the demineralizer;
- Be sure there is no water leakage ;
- Switch on the autoclave;
- Carry out one or several sterilization cycles to check for the correct connection and mainly to be sure of the absence of any leak.

At the end of every day always close the cock upstream the demineralizer

Connect the demineralizers exclusively to suitable autoclaves

NOTE: For the connection of the demineralizers to the autoclaves, see also the demineralizer manual.


FIG.29

14.3 Change of the filter

The autoclave indicates with a message on the display when the osmosis filter are getting exhausted. Replace the cartridge (*spare part 248-S-2*) as indicated in the demineralizer manual; After replacement, reset the cycles counter in the user setup (see chapter 8).

Sterilization means adopting precise operational protocol and applying a system. Autoclave is the final link of this system. Those who do not follow all of the various phases of the sterilization system will be unsatisfied with the final result. TECNO-GAZ S.p.A. produces all the devices for each phase of the sterilization cycle. Therefore it is necessary to:

1 - Prepare a precise written operational protocol to distribute to all operators assigned to sterilization.

TECNO-GAZ S.p.A. can provide these personalized protocols.

2 - Follow carefully and with attention all of the sterilization cycle phases:

DISINFECTION obligatory phase, in base of that is obligatory for the operator security, to be done with immersion in chemical liquids or thermo-disinfection.

CLEANING the most important phase which insures the removal of all types of residual, chemical and organic. The most suitable instruments are ultrasound baths.

DRYING essential phase to avoid corrosion of instruments and interference of the sterilization cycle.

PACKAGING essential phase for maintaining sterile the treated instruments.

STERILIZATION final phase of steam sterilization

16 INSTRUCTION FOR AUTOCLAVE SHIPPING

In case the autoclave must be shipped for repairs, controls, reactivation, revisions, validations follow the obligatory indications below:

- Use the original packaging; if this is no longer in your possession, use adequate packaging. The merchandise travels at risk to the sender.
- Ship the autoclave only (do not include any component contained in the accessories kit).
- Carefully clean the sterilization chamber and autoclave in general before shipping. In case it arrives dirty or with residual the autoclave will be returned without being repaired, or it will be put through a cleansing action and disinfection.
- Always empty the clean water tank through the attachment located on the back of the autoclave (*FIG.2-pos.2*)
- Always empty the used water tank through the attachment located on the back of the autoclave (*FIG.2-pos.5*)
- Indicate by letter and insert in the package a document which indicates precisely the irregularity or service desired.

ATTENZION:

All non original packaging which arrives will be disposed of.

Autoclave will be returned with new and original packing to insure maximum protection for your autoclave during shipping. Cost of packaging will be charged to client.

17 PROCEDURE FOR SERVICE AND ASSISTANCE

The repair service dept is an autonomous division of MEDILINE ITALIA SRL (TECNO-GAZ GROUP), comprising a complex and professional structure, in compliance with the procedures set in our *“Quality Assurance Manual”*.

- When the autoclave arrives at MEDILINE ITALIA SRL it is consigned directly to the service center, which codes and creates the technical assistance document.
- Within 48 hours, it is examined by technicians and technical analysis along with cost estimate is prepared for the service requested. If the autoclave is not delivered as requested in the chapter 16 instructions for shipping, above times can sensibly differ.
- Cost estimate will be sent to the distributor client, who will transmit it to the final client for vision and approval.
- After MEDILINE ITALIA SRL receives written acceptance of the cost estimate, the autoclave will be serviced and reshipped according to times indicated on the cost estimate.
- Ship at your expense, otherwise you will be billed for shipping.

IMPORTANT:

If after following the remedies, the problem persists, call telephone assistance at MEDILINE ITALIA srl +39 0522 - 94 29 97.

- 1) **DURATION:** the product is covered by a 12 month or 1200 cycles guarantee.
- 2) **EXCLUSION:** the guarantee excludes:
 - a) faults originating from non-performance of scheduled maintenance and due to negligence or improper use of the product by the User;
 - b) periodic checks and maintenance;
 - c) repair or replacement of parts subject to wear that are fragile or have an unpredictable life, unless they prove faulty at the time of their delivery;
 - d) faults caused by personnel or originating from technical personnel transfer and transportation;
 - e) faults or damage originating from improper use or operating errors;
 - f) faults or damage originating from contamination in the water and air infeed systems, extraordinary chemical agents or electrical events;
 - g) faults or damage originating from the use of detergents, disinfecting agents, sterilizing fluids or processes that are not described explicitly in the instruction manual for operation and maintenance;
 - h) natural color change of plastic components.
- 3) **INSTALLATION AND TESTING REPORT:** an essential prerequisite for guarantee validity is the return of the installation and testing report after all its sections have been filled in by the User and the Fitter. The report must be returned within 15 days after installation, otherwise the guarantee becomes null and void.
- 4) **LIMITS:** the guarantee entitles its holder to repairs or replacement of faulty components free of charge. The right to replace the entire machine is excluded. As for components either applied or added to TECNO-GAZ products or by third parties, having their own guarantee certificate, the conditions, limits and exclusions indicated in the relevant certificates apply.

Under no circumstances shall carrying out of one or more than one repair during the guarantee term change the guarantee expiry date.
- 5) **CLAIMS:** in the event of the Purchaser making a claim relating to the application of the guarantee and the quality or conditions of the delivered equipment, the said Purchaser shall in no way suspend and/or delay payments.
- 6) **DISPUTES:** if a dispute arises on the application and interpretation of this Certificate of guarantee, the Court of PARMA (Italy) shall be the court of jurisdiction, regardless of the location in which the purchase agreement was stipulated.
- 7) **DEROGATIONS:** any special derogations from the guarantee terms described in this section shall not imply the granting to the Purchaser of any right and shall be considered as granted for a specific case only.
- 8) **MISCELLANEOUS:** for any other issue, which is not covered by this Certificate of Guarantee, the rules described in the Civil Code and existing relevant Laws in the Republic of Italy apply.

The guarantee becomes null and void in the event that:

- a) the equipment presents damage caused by a fall, exposure to flames, liquid spillage, natural events, bad weather conditions or by any other cause that is not due to manufacturing defects;
- b) installation has not been performed in line with TECNO-GAZ instructions and by authorized personnel;
- c) the equipment is repaired, modified or tampered with by the Purchaser or by non-authorized third parties;
- d) when asking for an intervention under guarantee, the product serial number has been removed, deleted, counterfeited, etc.;
- e) the installation and testing report is not returned suitably filled in and signed within 15 days from installation;
- f) the Purchaser suspends and/or delays payment of amounts due, for any reason whatsoever, relating to the purchase and/or maintenance of the equipment;
- g) scheduled periodic maintenance or any other instruction envisaged in the manual for operation and maintenance is not carried out.

A

MAINTENANCE

Correct maintenance of the autoclave assures correct functioning of it and a secure saving in terms of time and costs for assistance and maintenance. The following operations are compulsory feasible by operators.

Cleaning of chamber

Every 20 cycles or once a week


Clean periodically the chamber, remove eventual deposits or debris, thus avoiding the introduction in the discharging circuit of material which can cause obstructions. For correct cleaning use only demineralized water and the abrasive sponge provided (*non abrasive side* - FIG.12).

To execute absolutely with cold chamber to avoid burns - Never use solvents, detergents, chemical solutions, descaling agents or other similar products.

Cleaning of chambre filter

Every 20 cycles or once a week

Pull up the filter (*spare code DXBA091*), paying attention not to damage it, wash with demineralized water and dry with cloth. Then replace in the seat, making sure that protrudes about 15mm.

Cleaning of trays and tray holder

Every 20 cycles or once a week

Clean with sponge supplied (non abrasive side) and demineralized water.

Change of the bacteriological filter

Every 200 cycles or when it assumes a dark color

Substitute the bacteriological filter (*FIG.2-pos.10*) turning it counterclockwise to unscrew and clockwise to screw. Use only original filter (spare code DAVA101).

Cleaning of the door gasket

Every 20 cycles or once a week

Periodically remove eventual residuals which are deposited on the circumference of the seal (*spare code DANA038*) using water and the sponge provided (non abrasive side), or else a moistened cloth.

ENGLISH

The following steps are recommended by the manufacturer. These operations must be performed only by trained technicians authorized by the manufacturer.

<i>Replacement of the water loading filter</i>	every 1 year / every 1000 cycles *
<i>Control of the water filter of flowmeter inlet</i>	every 1 year / every 1000 cycles *
<i>Cleaning Y brass filter downstream of the radiator</i>	every 1 year / every 1000 cycles *
<i>EV4 (water loading solenoid valve) cleaning</i>	every 1 year / every 1000 cycles *
<i>Vacuum pump valves replacement</i>	every 1 year / every 1000 cycles *
<i>Cleaning of the radiator</i>	every 1 year / every 1000 cycles *
<i>Fans cleaning</i>	every 1 year / every 1000 cycles *
<i>Control integrity of the pneumatic circuit</i>	every 1 year / every 1000 cycles *
<i>Control integrity of the electric circuit</i>	every 1 year / every 1000 cycles *
<i>Control thermal insulation of the chamber</i>	every 1 year / every 1000 cycles *
<i>Control tightening screws</i>	every 1 year / every 1000 cycles *
<i>Control temperature probes</i>	every 1 year / every 1000 cycles *
<i>Safety valve check</i>	Follow the guidelines of the manufacturer (see instructions)
<i>Control the state of wear of automatic lock</i>	every 1 year / every 1000 cycles *
<i>Level probes check</i>	every 1 year / every 1000 cycles *
<i>state of tank check</i>	every 1 year / every 1000 cycles *

* valid the condition of the reached first

Regional/national regulations fix the mandatory of autoclave validation that checks sterilization parameters.

Periodic validation must always be performed by trained technicians using appropriate measurement instruments.



In the case of non-mandatory, validation is a voluntary operation which aims to certify that the autoclave is fully functional and fit for sterilization of loads for which it was designed and built.

B SOFTWARE VISUALIZATION CYCLES


Insert the memory card on the front of autoclave into a memory card reader connected to your PC. Open "My Computer" (or similar) and double-click on external disk.

Open *index.htm*

Make sure that the ActiveX control is active to provide access to resources.

Press Apri file log and browse through the folders until you reach the one called LOG containing the save file of cycles made by the machine.

After selecting the desired cycle, in the mask are loaded the contents of all files present on the memory card containing the generic data cycles (date, cycle number, cycle type, outcome).

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Fine

Risorse del computer 100%

Selecting a cycle from a list, you will get a second mask containing all the details of the cycle.

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
ENGLISH

Press Anteprima di stampa, will be displayed a page in text format showing autoclave data and autoclave data cycle in a structure suitable to be printed.

Riga	Data	Time left	T1	T2	P1	V rete	Cycle step
1	17/12/2009 14:37:48	0.00	73.74	75.34	-0.011	234.7	START-UP
2	17/12/2009 14:37:56	0.00	73.97	75.47	-0.010	232.8	START-UP
3	17/12/2009 14:37:56	960.00	74.47	75.86	0.000	234.6	VACUUM
4	17/12/2009 14:38:05	960.00	74.55	76.15	-0.151	226.5	VACUUM
5	17/12/2009 14:38:15	960.00	74.76	76.85	-0.210	227.4	VACUUM
6	17/12/2009 14:38:26	960.00	75.70	77.56	-0.342	223.5	VACUUM
7	17/12/2009 14:38:36	960.00	76.09	77.74	-0.400	225.8	VACUUM
8	17/12/2009 14:38:46	960.00	75.57	77.09	-0.483	225.1	VACUUM
9	17/12/2009 14:38:56	960.00	74.37	75.94	-0.542	226.4	VACUUM
10	17/12/2009 14:39:05	960.00	72.87	74.00	-0.608	225.8	VACUUM
11	17/12/2009 14:39:15	960.00	72.25	73.06	-0.638	226.7	VACUUM
12	17/12/2009 14:39:26	960.00	69.47	70.00	-0.690	225.9	VACUUM
13	17/12/2009 14:39:36	960.00	66.98	67.40	-0.718	226.7	VACUUM
14	17/12/2009 14:39:46	960.00	65.91	66.30	-0.728	226.7	VACUUM
15	17/12/2009 14:39:56	960.00	64.92	65.31	-0.739	224.0	VACUUM
16	17/12/2009 14:40:06	960.00	64.24	64.60	-0.744	227.8	VACUUM
17	17/12/2009 14:40:17	960.00	63.95	64.14	-0.760	224.9	VACUUM
18	17/12/2009 14:40:27	960.00	62.86	62.44	-0.777	225.4	VACUUM
19	17/12/2009 14:40:37	960.00	62.63	61.97	-0.776	226.8	VACUUM
20	17/12/2009 14:40:47	960.00	62.29	60.93	-0.793	226.7	VACUUM

Selecting: File->Stampa you can select the desired printer and start the process of printing of the document.

Before removing the memory card you should disconnect the device, to avoid errors or file system corruption on the card itself.

We recommed to make a copy of the entire contents of the memory card on your computer at least once a week and arrange for appropriate storage backups

How to read the press of a cycle

Matricola: AUR00X090014

Serie firmware: T1A0000

Rel. firmware: 2

N. ciclo: 10

Desc. ciclo: Ciclo 121°C

F
G

Riga	Data	Time left	T1	T2	P1	V rete	Cycle step
1	17/12/2009 15:22:11	0.0	47.79	47.79	0.1	213.3	START-UP
2	17/12/2009 15:22:18	0.0	47.84	47.84	0.5	213.9	START-UP
3	17/12/2009 15:22:18	1800.00	47.84	47.84	0.5	213.9	VACUUM
4	17/12/2009 15:22:25	1800.00	47.72	47.72	-0.5	213.9	VACUUM
5	17/12/2009 15:22:56	1800.00	47.90	47.90	-0.1	206.0	VACUUM
6	17/12/2009 15:23:25	1800.00	48.71	48.71	-0.5	206.9	VACUUM
7	17/12/2009 15:23:56	1800.00	49.71	49.71	-0.3	208.3	VACUUM
8	17/12/2009 15:24:26	1800.00	51.09	51.09	-0.7	208.1	VACUUM
9	17/12/2009 15:24:49	1800.00	52.13	52.13	-0.9	209.0	VACUUM
10	17/12/2009 15:24:52	1800.00	52.18	52.18	-0.784	214.6	VACUUM
11	17/12/2009 15:25:00	1800.00	60.31	60.31	-0.784	212.7	HEATING

A	Indication of date and time when the machine makes a reading of parameters
B	Time remaining at the end of cycle (in seconds – is the sum of the time of sterilization and the drying)
C	Reading the two temperature probes
D	Reading of chamber pressure
E	Reading of voltage
F	Cycle information, serial machine and firmware version
G	Operational phase of the cycle

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EN

This manual must always be kept with the product, in complying with the Directives of European Community. TECNO-GAZ reserves the right to modify the enclosed document without notice. TECNO-GAZ reserves the property of the document and forbids others to use it or spread it without its approval.